


The Critical Role of Bench Tests in Automotive Lubricant Specifications, Development and Production Control

*How Tannas & King Instruments
Serve the Industry*





Midland, Michigan USA

**Independent
Testing
Laboratory**

1969



SavantLab.com

**Laboratory Testing
Instruments**

1981



TannasKing.com

1994



**Global
Engine Oil
Database**

1984



InstituteOfMaterials.com

**Precision
Machining**

2019



ExcellMfgInc.com





**Midland,
Michigan**





Founder & Inventor
(1928 – 2022)

**Theodore (Ted)
Selby**

(58 Patents)



50th Anniversary Celebration

1969 - 2019



KING
DESIGNATION

Tannas & King Philosophy

Serve
CUSTOMER
Needs

- ◆ Produce the finest available equipment
- ◆ Solve critical problem areas in industry
- ◆ Training & consultation

Serve
INDUSTRY
Needs

- ◆ Bench test development
- ◆ Critical information production
- ◆ High correlation with engine tests & field performance



Instrument Innovations - Tannas




Quantum® Oxidation Tester

- Non-liquid dry cylinder
- RPVOT (D2272, D2112, IP 229, SH/T 0193)
- TFOUT (D4742)
- Grease Oxidation testing (D942)



TBS HTS Viscometer 2100E-F / 3000

- Tapered Bearing Simulator
- 40°C to 150°C
- 40+ position Auto-Sampler
- ASTM D4683, D6616, CEC L-036, SH/T 0618



TFAB Tannas Foam Air Bath

- Non-liquid dry bath
- Direct Drive motor - quiet
- Touchscreen interface
- ASTM D892, D6082, IP 146, GB/T 12579, SH/T 0722




Direct Cool II (non-liquid)

- Low temperature pumpability
- Gelation Index
- 90°C to -40°C
- ASTM D5133, D7110, SH/T 0732
- ASTM D2983, Procedure D (Automated)



KING
REGENERATION

Instrument Innovations - Tannas



**Tannas
Noack S2[®]**
(non-Wood's Metal)

- Evaporation Loss
- 'Tunable'; Calibrate to lab environment
- Collection of volatiles – PEI
- ASTM D5800-Procedure **D**, SH/T 0059



TEOST[®]
**Thermo-oxidation
Engine Oil Simulation Test**

- High Temperature Deposit Control
- TEOST[®] 33C → ASTM D6335, SH/T 0750
- TEOST MHT[®] → ASTM D7097, SH/T 0834
- TEOST Turbo[®] → ASTM D8447



SBT[®] +2
*(Multi-Purpose Low-Temp
Liquid Bath)*

- Easily replaceable Insert Modules for each test
- ASTM D5133, D2983, D7110, D4684, D445, D97, D2500, SH/T 0732



Instrument Innovations - King



BLB

Brookfield Liquid Bath

- SimAir™ Test Cells
- Independent Sample Analysis
- 30°C to -70°C
- ASTM D2983, DIN 51398



MRV TP-1

Mini-Rotary Viscometer

- Low-Temp Pumpability & Yield Stress
- Direct Refrigeration Technology
- ASTM D3829, **D4684**, D6821, D6896, SH/T 0562





Tannas & King Refrigeration

Lab Instruments

50+ Distributors

Worldwide

Instruments sold to over

70 countries

Instrument Sales Worldwide

2022



Industry Specifications

**SAE
J300**

**Engine Oil Viscosity
Classification**

API & SAE Viscosity Grade
Symbol



ILSAC

**International Lubricant
Standardization &
Approval Committee**

GF-2 thru GF-5+, GF-6(a&b)
API S? Category equivalent



OEM

OEM-Specifications

ACEA -Europe



IFC



General Motors



Industry Specifications

ATF

**Automatic Transmission
Fluid:**

**GM Dexron®
Ford Mercon®
Chrysler (FCA) ATF+4
OEM Genuine**

Industrial Oils

**Wind & Gas Turbine
Manufacturers:**

**GE
Siemens
Vestas
Goldwind
United Power
Kawasaki Heavy
Alsald Energia**

China National Standard

**GB 11121
Gasoline Engine Oils**

**GB 11122
Diesel Engine Oils**

**New D1 Specification for
Diesel**


**New Gasoline & EV
Specification Coming**



SAE J300 Viscosity Classification



**Tannas/King
Instruments
included in
J300**



SAE Viscosity Grade	Low Temp. (°C) Viscosity, cP		High Temperature (°C) Viscosity		
	Cranking, Max. (CCS)	Pumping Max. (NYS) (MRV TP-1)	Kinematic (cSt) @ 100°C Min.	Kinematic @ 100°C Max.	High-Shear (cP) @ 150°C D4683, D4741 and D5481
0W	6200 @ -35	60,000 @ -40	3.8	--	--
5W	6600 @ -30	60,000 @ -35	3.8	--	--
10W	7000 @ -25	60,000 @ -30	4.1	--	--
15W	7000 @ -20	60,000 @ -25	5.6	--	--
20W	9500 @ -15	60,000 @ -20	5.6	--	--
25W	13,000 @ -10	60,000 @ -15	9.3	--	--
8	--	--	4.0	< 6.1	1.7
12	--	--	5.0	< 7.1	2.0
16	--	--	6.1	< 8.2	2.3
20	--	--	5.6	< 9.3	2.6
30	--	--	9.3	< 12.5	2.9
40	--	--	12.5	< 16.3	3.5 (0W-40, 5W 40, 10W-40)
40	--	--	12.5	< 6.3	3.7 (15W-40, 20W-40, 25W-40, 40)
50	--	--	16.3	< 21.9	3.7
60	--	--	21.9	< 26.1	3.7

ILSAC GF-6 (A & B) Engine Oil Specifications


SAE J300 Multi-grade Oils:

- CCS (D5293)
- **MRV** (D4684)
- KV @ 100°C (D445)
- **HTHS @ 150°C** (D4683, D4741, D5481)
- **SBT®/ Gelation Index** (D5133)

Engine Test Requirements:

- Oil Thickening (Sequence IIH)
- Sludge and Varnish Test (Sequence VH)
- Valvetrain Wear (Sequence IVB)
- Bearing Corrosion (Sequence VIII)
- Fuel Efficiency (Sequence VIE)
- Chain Wear, Low Speed
- Pre-ignition Prevention

ILSAC
Starburst & Shield




Bench Test Requirements:

- Phosphorous % (D4951)
- Phosphorous Volatility (D8111)
- Sulfur Content (D4951, D2622)
- **Volatility (D5800 B/D)**
- **TEOST 33C (D6335)**
- Filterability (D6794, D6795)
- **Foam Tendency (D892, D6082 high temperature)**
- ROBO (D7528) or Sequence IIHA Test (D8111)
- Shear Stability (D6709)
- Homogeneity (D6922)
- Ball Rust Test (D6557)
- Emulsion Retention (D7563)
- Elastomer Compatibility



OEM Engine Oil Specifications

Tannas/King instruments included in:
ILSAC, ACEA, IFC, dexos



European
Automobile
Manufacturers
Association





Tapered Bearing Simulator (TBS) Viscometer

HTHS Viscometer
TBS 2100E-F

*The instrument that opened the door to
understanding engine dependence on
high shear viscometry*



High Shear Rate Viscometry - Automotive Oils

Newtonian Oil

- Straight weight
- SAE 30
- No Viscosity Modifier
- Viscosity = Shear Stress/Shear Rate

Non-Newtonian Oil

- Multi-grade
- 5W30
- Viscosity Modifier (VM)
- Viscosity decreased with increasing shear rate (temporarily)



High Shear Rate Viscometry


Multigrade oils benefits, easier engine cold starting with acceptable viscous protection of operating engine.

As **multigrade** engine oils increased, it was found they had interesting viscometric properties that required **new understanding.**

New area of viscous measurement was developed to further this understanding – **HTHS** testing.



High Shear Rate Viscometry



High shear rate viscosity influences whether sufficient hydrodynamic lubrication occurs to protect opposing moving surfaces.

Viscosity

is a form of
friction

A resistance to flow caused by the internal molecular interaction of the oil.

Viscosity **protects** the engine from wear,
but is a major source of
fuel consumption.

High Shear Rate Viscometry

Viscosity Losses

are caused by viscosity contributions of **VI Improvers** before and after shear degradation.

High Temperature, High Shear

rate viscosity predicts how well the oil will protect engine bearings at very high shear rates.

The **TBS Viscometer** has generated the most published information.



History of TBS HTHS Viscometer

In 1979, Exxon's Paramins division commissioned the development of a very high shear rate, relatively high temperature (HTHS) viscometer with commercial development rights to **Savant Inc.**

Tannas Co. was formed in 1981 to manufacture and market **first commercially** available HTHS device... the rotational **TBS** Viscometer.

The **importance** of TBS became evident with ASTM studies correlating engine bearing oil-film thickness (BOFT) with high shear viscometry.



TBS Viscometer Basics



Tapered Rotor & Stator Geometry

- High speed rotation
- 3.5 micron gap (~1/30 human hair)
- Change shear rate (rotor position & speed) during operation




Electric Thin-Film Heater

- Constant temperature control
- 40°C to 200°C $\pm 0.1^\circ\text{C}$
- Rapid temperature recovery between measurements


Precise Temperature Control



TBS Rotor & Stator




High Shear Testing Fundamentals



Rotational HTHS Viscometer

simulates engine bearing lubrication conditions



Capillary Viscometer

--does not simulate engine
--adds flow complexity with parabolic flow behavior



TBS 2100 $E-F$


12 Speed Selections: **800-8,000 RPM**

Shear rates: **50,000 to $7 \times 10^6 \text{ sec}^{-1}$**

Chase Flush Sample Injection

42-position AutoSampler with PC controlled calibrations, injections, etc.

Only **'Absolute'** HTHS Viscometer



NEW

TBS 3000 HTHS Viscometer

*Customer feedback and
advancing technology
lead to new redesigned
TBS model for
80°C, 100°C, and 150°C
operation.*



High Shear Rate Viscometry

**Dual-ring 40-position
Autosampler**

**Reduced sample
volume per test
(40 mL)**


**10 easily replaceable
Cal. Oil Containers**

**Remote Access via Ethernet/
Connection to LIMS**



High Shear Rate Viscometry

HMI Touchscreen with TBS Touch Integrated Software.



High Shear Rate Viscometry




**Integrated Custom
Temperature Control
System (80°C to 150°C)**

**No computer, DAQ
Board, reducing
Cabling**

Auto Rubbing Contact

**Custom LED lighting
based on operating
conditions**

TBS Viscometer Models



Semi-Auto TBS 2100E-F
Full-Auto w/Manual Injection



Full-Automation TBS 2100E-F



New TBS 3000



TBS Industry Status

ASTM D4683

- 150°C at 1×10^6 reciprocal seconds
- JPI-5S-36-03 (Japanese)

ASTM D6616

- 100°C at 1×10^6 reciprocal seconds
- Required for railroad engine oils & dexos™ specs.

CEC L-036

- European method
- 100°C and 150°C operation at 1×10^6 reciprocal seconds
- SH/T 0618

SAE J300

- Referee Method = D4683, 150°C at 1×10^6 sec⁻¹
- ILSAC GF-1 to GF-6

GB-11121 GB-11122

- Chinese National Standards (Gas & Diesel)
- via SAE J300 & SH/T 0618


Advanced TBS Applications

Multi-Temperature Analysis




Multi-Shear Rate Analysis


**Viscosity Loss Trapezoid
(Profile) - VLT**




Multiple Temperature Viscometry



Smooth viscosity-temperature data over broad temperature range




Multiple Shear Rates



Broad shear rate range:

- Evaluates properties and types of VI Improvers based on their response to very high shear rates





Comments or Questions?

explore our new website
TannasKing.com



Oxidation Stability of Lubricating Oils

Quantum® Oxidation Tester
RPVOT
TFOUT
Grease D942
Advanced Applications




Introduction – RPVOT

(Rotating Pressure Vessel Oxidation Test)

Originally designed to simply compare in-line samples of turbine oils.

Later the test developed into testing fresh and in-service turbine and engine oils.

This test eventually developed limits and is in several specifications.



Images: www.Pixabay.com

Introduction – RPVOT

(Rotating Pressure Vessel Oxidation Test)

The RPVOT per ASTM D2272 / IP229 is the primary tool used in the analysis and prediction of turbine oil life.


Parameters

- Sample Size: 50g
- Temperature: 150C
- Initial O₂ Pressure: 90 psi
- Pressure increases to ~190 psi
- Copper Catalyst Coil
- 10 mL of water (oil sample)
- Rotational Speed: 100 ±5 rpm
- Rotational Angle: 30°
- Test monitored by Computer or Strip Chart
- Pressure drop of 25.4 psi signals end of test (Break-point).



Introduction – RPVOT

(Rotating Pressure Vessel Oxidation Test)



Introduction - RPVOT

RPVOT Liquid baths utilize a high-temperature oil bath (often silicone)

Silicone vapors are a health hazard and can be very expensive to replace

RPVOT liquid baths use multiple pressure vessels (requiring considerable hood space).

RPVOT liquid baths → messy, smelly, & time-consuming to replace & repair



Quantum[®] RPVOT Oxidation Tester

Developed in response to customer needs for safe and efficient means of running RPVOT tests.

Single position 'Dry Cylinder' heated stainless steel oxidation vessel with pressure transducer.

Convenient front-loading of samples, oxygen charge & release valves. *(No Hood Requirement)*




Quantum[®] RPVOT Oxidation Tester

Designed for either multiple (side-by-side) bench-top alignment or "stand-alone" operation



Automation Package capable of monitoring and recording up to four (4) *Quantum[®] units independently.*



Grease Oxidation Testing (D942)



Grease Oxidation test D942 is used with a special tipping Rack.

This Rack allows the unit to sit at 90° in order to load the Grease Tray.



Quantum[®] Industry Status

ASTM D2272

- IP229; SH/T 0193
- RPVOT -- Oxidation stability of Turbine Oils

ASTM D2112

- Oxygen Stability of inhibited mineral insulating oils
- 140°C @ 90 PSI oxygen & all other D2272 parameters

ASTM D4742

- TFOUT -- Oxygen Stability of engine oils
- 160°C @ 90 PSI oxygen
- TFOUT Catalyst Package (Metal & Fuel Catalyst, Water)

ASTM D942

- Oxygen Stability of Grease
- 99°C @ 110 PSI oxygen
- 20g of grease in 5 glass dishes (4g per dish)



Advanced *Quantum*® Applications



**Sample
Temperature
Monitoring**



**Sample Extraction
& Injection**



Sample Extraction/Injection Device



New

Available NOW!

Quantum[®] Pro
Oxidation Stability of
Lubricating Oils



New Quantum[®] Pro

Automated Fill and Purge Control

Integrated Water Calibration


Live Graph of Test Data

Integrated Sample Temperature Probe

Lower Power Usage (20 - 30%)

USB Ports to Transfer Records to Flash Drive or PC

Automated Test Report Submission Over Ethernet





Comments or Questions?

explore our new website
TannasKing.com





Foam Tendency & Stability of Lubricating Oils

Tannas Foam *Air* Bath
TFAB




Introduction

Because foam is mostly entrained air, and air is a poor heat conductor, foam interferes with the fluids' cooling performance. Moreover, air is a poor lubricant.

Tendency of oils to foam can cause inadequate lubrication, cavitation, and overflow loss of lubricant, resulting in mechanical failure.

Antifoaming agents are used to reduce and control the amount of foam, and this foam test is a gauge to determine the effectiveness of the additive.




Tannas Foam *Air* Bath

Bench-top AIR BATH – safe & clean alternative to messy liquid baths

DIRECT DRIVE motor for quiet, maintenance-free operation



- 6-sample load capacity
- Touchscreen controller

Tight Temperature Control
Cooling Connections to water supply




TFAB Overview

Touchscreen
Controller



TFAB Parts





Diffuser Stones w/Centering Ring



Diffuser Stone with Certificate

Recommended TFAB Testing Order



TFAB Industry Status

ASTM D892 / IP146

- Foaming Characteristics of Lubricating Oils
- Sequence I, II, III (24C & 93.5C)

ASTM D6082

- High Temperature Foaming Characteristics of Lubricating Oils
- Sequence IV (150C)

ASTM D1881 & D7840


- Foaming Tendencies of Engine Coolants
- Foaming Tendencies of Non-Aqueous Engine Coolants

Chinese Standard

- GB/T 12579 = D892
- SH/T 0722 = D6082
- SH/T 0066 = D1881

ILSAC GF specs

- Many other industry specifications





Comments or Questions?

explore our new website
TannasKing.com






Engine Oil Volatility... w/Expanded Capabilities

Tannas Noack S2®



Noack Volatility Test




Noack History

Oil formulators have long-since recognized the value of determining the volatility of automotive lubricants.

However, the only test considered sufficiently meaningful was the Noack test that had been developed by Dr. K. Noack in the 1930s.

This test used Wood's Metal as a heat-transfer medium – a toxic substance which had carcinogenic properties – a seriously negative factor.

Despite this, in the mid-1990s it became an ASTM Test Method and thus a standard in choosing base stocks and finished lubricants.





History

During this timeframe, Savant Laboratories in cooperation with Tannas Co. developed an alternative Noack volatility approach.

A new instrument that eliminated the toxic Wood's Metal heating and also collected the volatiles was conceived and patented.

In addition to gathering Noack volatility data needed for industry specs., this innovation permitted study of the volatilized oil and comparison to the original and residual oil (with some interesting and important findings).

The original Selby-Noack™ is ASTM D5800C and the new Tannas Noack S2® is included as D5800D (which matches Procedure B results) and CEC L-40.




Noack Volatility Test

ASTM D5800

PARAMETER	VALUE
Sample Size	65 grams
Operating Temperature	250°C
Vacuum Control	20 mm of H ₂ O
Test Time	1 Hour
Result	% weight difference



Noack D5800 Volatility Test



**Tannas Noack S2[®] → Included as D5800
Procedure D
(matches Procedure B results)**



Tannas Noack S2[®] Volatility Test


Redesigned for easier operation & improved precision

- ➔ Metal Reaction Vessel
- ➔ Quick Connect Fittings

Interchangeable Orifice Tubes to 'Tune'/calibrate to lab environment

Non-Wood's Metal with same D5800 test parameters

Collects volatilized materials for PEI



Volatilized Material Analysis

This instrumental approach to the Noack volatility test permitted collection of nearly all volatilized material. Used to calculate the Phosphorus Emission Index (PEI). Characterizes phosphorus emission potential into exhaust stream.



Clean Coalescing Filter



After Test



Tannas Noack S2[®] – Parts

Reaction Vessel Assembly



Tannas Noack Industry Status

ASTM D5800

- Procedures A & C correlate (Annex). Procedure B & D correlate – and to Procedures A and C with correction factors.
- SH/T 0059 Chinese Standard

API / PC-11

- SJ | SL | SM | SN | SP Categories
- 15% Spec limit

ILSAC

- GF-2, GF-3, GF-4, GF-5+, GF-6 (A&B) Specifications
- 15% Spec limit (Proc. B or D only)

CEC L-040

- Fully accepted and included in early 2022.

dexos™

- Gen 1 & 2 = 13.0% limit (avg. of 3 runs)
- Gen 3 = 12.5% limit (avg. of 3 runs)



Expanded Capabilities

Historically, 250°C was intended to expedite the test. Recently, industry members are realizing this temperature may be too high with questionable correlation to engine oil consumption for modern engines and EV powertrains.



Industry partners have been working with Savant Labs to test at a variety of temperatures and times.



Recent data has shown that 150°C is the average temperature for critical engine areas.

OEMs are considering Noack testing at 150°C @ 12 hr. run time for future specifications.





Comments or Questions?

explore our new website
TannasKing.com



Low Temperature Behavior of Engine Oils

Pumpability

Gelation Tendencies



ASTM D5133
SH/T 0732



ASTM D4684
SH/T 0562



Low Temperature Behavior

Low-temperature **STARTABILITY**

Modern Engine designs,
fuel injection

is a problem of the
PAST

modern
**HYDRO-TREATED
BASE OILS**



the need for special
**POUR POINT
DEPRESSANTS**



**VI
IMPROVERS**



Controlling **PUMPABILITY** at
LOW SHEAR STRESS is now a
challenging area of
oil & additive development



Low Temperature Pumpability Failures

Flow-limited behavior

- ◆ Too viscous to pump (milk-shake w/straw)
- ◆ An inconvenience
- ◆ Viscosity reduced by increasing temperatures

Air-binding behavior

- ◆ Combination of gelation and viscous flow
- ◆ Structure (component nucleation/crystallization)
- ◆ Can occur at any low-temperature
- ◆ Catastrophic to operating engine



Air-Binding Behavior

Low-Temperature Engine Tests

'Air-Binding' Process

Oil supply lines to
valve-train and crankshaft


Air entrained in lines

Atmospheric
Pressure

Oil pump

Gelated oil
in sump

Oil Screen



When oil
GELATES
in the sump,

the **OIL PUMP**
draws a **core** of
structured oil

forming an
AIR TUBE

The **OIL PUMP**
becomes

AIR-BOUND

causing
**ENGINE
FAILURE**

Low Temperature Pumpability

ASTM
D3829

1979

- ◆ Mini-Rotary Viscometer (16hr cooling)
- ◆ Correlated to cold-room engine tests (BPT)
- ◆ Thought D3829 could predict low-temp pumpability of engine oils

1980-81

- ◆ Massive field failures in Europe & North America
- ◆ D3829 not correlative to natural cooling cycles



Low Temperature Pumpability

- ◆ Imitated field conditions
- ◆ Exact cooling condition to pre-set the gelated oil structure

**Special
Cold-Room
Engine
Tests**

- ◆ Became ASTM D4684
- ◆ SAE J300 | eventually ILSAC GF series
- ◆ TP-1 conditions the sample using a variable programmed cool rate (45+hrs)

**MRV
TP-1**



King MRV TP-1

Direct-refrigeration system (no circulating bath) (Procedure A of D4684)

Tight temperature control ($\pm 0.1^{\circ}\text{C}$)

(10) Sample Capacity with Removable test cells

New rotor design & cooling air flow
(limits water & ice interference)

**ASTM D4684 & D6821 (Driveline)
and D6896 (used Diesel Oils); SH/T 0562**



Scanning Brookfield Technique (D5133)

1982

SCANNING BROOKFIELD TECHNIQUE (SBT®)

- ◆ First pumpability bench test to correlate to field failures
- ◆ Forced oils to gelate (if they would)

**The SBT®
Provided:**

**Constant
laminar
motion**
(0.3 rpm)

**Precise,
gradual
cooling**

**Continuous
viscosity
measurement**

**Revealed onset
and degree of oil
Structure-
building**



ASTM D5133


Slow cooling profile → 1°C/hour (-5 to -40°C)

Precision made titanium rotor & glass stator

Specially modified digital viscometer

Constant recording of sample viscosity

Scanning Brookfield Rotor/Stator Cell




SBT[®]

ASTM D5133


Automation program

Real-time display of data

Program prints composite plots of results



SBT[®] (D5133 & D7110)




Flow-Limited
viscosity related
behavior



Exposes more
subtle
air-binding
tendencies
of field-failing oils.

SBT[®] (D5133 & D7110)



Obtain the
Gelation Index*

to measure the degree of gelation.

* The first derivative of the MacCoull, Walther, Wright equation



Only technique that
generates the
GELATION INDEX
and
GELATION INDEX
TEMPERATURE

First procedure to show
field failures
were
OIL-RELATED
rather than
engine-related

SBT
correctly predicted
ALL
field-failing oils

SBT
shows high correlation
on both
AIR-BINDING &
FLOW-LIMITED
oils from ASTM
cold-room tests

SBT[®]

D5133

D7110

Most complete
analysis available on
engine oil
pumpability.

SBT[®] & MRV



Engine Oil
specifications
require both
SBT & MRV TP-1
methods



KING
REFRIGERATION

(2) SBT[®] Instrument Models



Direct Cool II



SBT+2 Multi-purpose Bath



Direct Cool II

SBT®
(D5133, D7110)

Low-Temp Viscosity
(D2983, Procedure D)

*Customer feedback and
advancing technology
lead to new redesigned
Direct Cool II model.*



Direct Cool II

Thermal Electric Cooling (efficient, more robust, no refrigerants)

All Digital Output

Automation Package (up to 4 units)

Touch-screen Controller (easy navigation of test parameters)

Self-Contained Pre-Heating to 90C

Eliminate DAQ Board, Universal Power

Multi-Purpose: Automated D2983, Procedure D



SBT+2 (multi-purpose bath)

Multi-Purpose Low-Temp Bath

Programmable & Manual Control

Cooling to -70°C

Insert Modules for each Test

ASTM D5133, D7110, D2983, D445, D97, D2500



SBT+2 – Other Test Insert Modules

**D2500, D97
Cloud & Pour Point**



**D2983-Procedure C
Brookfield Viscosity**



**D4684
Mini-Rotary Viscometer**

**D445
Kinematic Viscosity**



SBT[®] Industry Status

ASTM D5133
ASTM D7110

- D5133 → cooling profile: 1°C/hour from -5°C to -40°C
- D7110 → cooling profile: 3°C/hour from -5°C to -40°C

API

- SL
 - SM
 - SN
- SP** } Engine Oil Categories

dexos[™]

- General Motors specification

ILSAC

- GF series Passenger Car Specifications

GB-11121


- Chinese National Standard
- SH/T 0732



Low Temperature Behavior of Engine Oils

Low-Temperature Viscosity Testing (formerly known as Brookfield Viscosity)

ASTM D2983, Procedure C
IP 267
DIN 51398



Low Temperature ATF Pumpability

Low-Temperature Viscosity





Image: © ASTM D2983-15, pg. 4

1950's


- Automatic transmission failures in cold temperatures
- Auto manufacturers began low-temperature studies

1971

- Procedure written and accepted as ASTM D2983
- Balsa wood block & cold Air cabinet
- Samples removed & tested with a rotational, bench-top Viscometer (Brookfield originally used)



Liquid Baths with Tannas SimAir® Test Cell




A modified stator incorporates:

- A thin, gas-filled Dewar section
- This modulates heat transfer (cooling rate) from test fluid to liquid bath

1990's

- Tannas developed and patented the SimAir® Test Cell
- Used in a constant temperature liquid bath

- The SimAir® Test Cell led to dramatic improvements with:
 - Test Operation
 - Usability
 - Precision



SimAir[®] Test Cell


Innovation for ASTM D2983

No duplicate sample runs required.

Simply set the liquid bath at the desired final temperature.

Insert sample for 16-hour exposure.

Dewar stator cell assembly



BLB Innovation

SimAir[®] Test Cells


Provides proper sample temperature control

Used exclusively in King BLB or Tannas SBT+2 liquid baths

Each Test Cell has its own independent cooling profile


Each cell can be added or removed from the bath at any time

Increased productivity and accuracy





Sample Temperature, °C



Comparison of a single wall stator and #4 spindle to the SimAir® Cell.

Note the effectiveness of a constant temperature liquid bath in simulating air bath results when using the SimAir® Test Cell.

Low-Temperature Viscosity *(Air baths – Liquid baths – SimAir®)*

Air

Transfer samples from air bath to Benchtop Viscometer

◆
Balsa wood carrier

Air to Liquid

Transfer samples from air bath to constant temperature liquid bath 30 minutes before analysis

◆
Viscometer mounted on bath

Programmable Liquid Bath

Energetic Refrigeration System

◆
Attempts to emulate cooling rate of air bath

◆
Force cool each sample the same

SimAir® Test Cell

Simpler approach

◆
Constant temperature liquid bath

◆
Simulates cooling rate of samples as in air bath

◆
Each Cell independent of others in bath

Liquid Baths with Tannas SimAir® Test Cell

King BLB Liquid Baths


Holds up to 12 samples
Broad Temperature Range
($\pm 0.1^{\circ}\text{C}$ control)

BLB 701: $+30^{\circ}\text{C}$ to -40°C

BLB 702: $+30^{\circ}\text{C}$ to -70°C



King BLB-DIN Liquid Bath



Meets ASTM D2983 & DIN 51398

Constant temperature for ASTM (SimAir®)

Programmed for DIN (1°C/min) (Single wall stator)



+30°C to -55°C

Holds up to (12) samples

SBT+2 with D2983 Insert

SimAir[®] Test Cells

Holds up to (8) samples



KING
INFORMATION



Comments or Questions?

explore our new website
TannasKing.com



Improving Prediction of Deposits

High Temperature Deposits

TEOST® 33C

TEOST MHT®

TEOST® Dual

TEOST Turbo® /
(3 tests, 4 models)

NEW



TEOST[®] 33C for *Turbocharger Deposits*

1985

Savant understood
**DEPOSIT
CONDITIONS**
fall into two
categories →

- 1: Precursors at high-normal engine temps
- 2: High temp deposit-inducing zone

1989

Cooperative
efforts with
CHRYSLER
Corporation


Development
of bench test
to predict the
**DEPOSIT
PROTECTION**
of an oil



TEOST® 33C for *Turbocharger Deposits*

- Exhaust heat from turbine shaft (600-700°C) is absorbed by the engine oil pumped through the bearings.
- Earliest stages of coking occur just after the engine is shut off.
- Heat damages the oil when engine oil circulation is stopped causing deposits to form, blocking oil passage and leading to failure.

Cut-away Schematic of Turbocharger



TEOST[®] 33C for *Turbocharger Deposits*


2-hour performance test

Steel depositor rod heated cyclically between 200°C & 480°C for 12 cycles


Rod weighed before & after test to determine deposit levels

Deposit level of rod + filtered oil (mg)

- 60 mg for GF-2 and 30 mg for GF-5 & GF-6



TEOST[®] 33C Cyclic Heating Conditions



Based on oil temperatures in the turbocharger shaft during operation and cool-down.

Depositor Rod
Total Heat Cycle
Duration:

**9
minutes**

Depositor Rod
Heat Cycle
Range:


**200°C
to
480°C**

Number of
Heat Cycles
per Test:

12

33C Operating Parameters

PARAMETER	VALUE
Total Oil Volume	116 mL
Oil Pumping Rate	0.40 g/min
Reactor Temperature	100°C
Moist Air Flow	3.5 mL/min
N ₂ O Flow	3.5 mL/min
Iron Naphthenate	193 microliters
Position of controlling Thermocouple in Rod	~70 – 75 mm from top of Rod



TEOST® 33C Industry Status

ASTM D6335

- Part of ASTM Test Monitoring Center Program

ILSAC

- **GF-2, GF-5+ & GF-6** Engine Oil Specifications
- GF-2 = 60 mg; GF-5+, GF-6 = 30 mg

API

- **SJ, SN & SP** categories for engine oils
- ASTM D4485

GB-11121

- SH/T 0750

Development of TEOST MHT®

TEOST MHT® bench test developed for Chrysler to predict **PISTON VARNISH** & **DEPOSIT TENDENCIES** of engine oils.

MHT developed to correlate with **Peugeot TU3MH**
(a European piston varnish & deposit rating test)

Correlation using the MHT was considered good, with:

- ◆ **Specially modified pistons** (required by engine test) **AND**
- ◆ **Commercially available pistons**

TEOST MHT[®] for *Engine Deposits*

24-hour performance test

Wire-wound depositor rod (Specially treated)

- Rod held at 285°C
- 8.5 g of oil recirculating at 0.25g/minute


MHT Catalyst – represents concentration of certain metals in used oil

Rod & Filter weighed before & after test to determine deposit levels



TEOST MHT[®]

Depositor Assembly For the TEOST MHT-4 Protocol



Oil flows down spiral path of wire


Air flow of 10 mL/minute

Volatile oil removed & collected

Deposit level of rod + filtered oil (mg)

MHT[®] Operating Parameters

PARAMETER	VALUE
Total Oil Volume	~ 8.5 g
Oil Pumping Rate	0.25 g/min
Operating Temperature	285°C
Test Time	24 Hours
Dry Air Flow Rate	10 mL/min \pm 0.2
MHT [®] Catalyst	~0.0131 g/g of Oil



TEOST MHT® Industry Status

ASTM D7097

- Part of ASTM Test Monitoring Center Program

ILSAC

- **GF-3, GF-4 & GF-5+** Engine Oil Specifications
- 35 mg limit spec in **GF-5+**

API

- **SL, SM & SN** categories for engine oils
- **ASTM D4485**

GB-11121

- **SH/T 0834**



Updating the Turbocharger Bench Test

To align with new turbocharger technology, development of an updated bench test began in early 2017 for key reasons:

- **To assist with development of new oil formulations**
- **For screening engine oil formulations prior to running updated turbo engine tests**
- **To provide on-going evidence of engine oil quality**



TEOST Turbo[®]: *Experimentation & Design*

The TEOST Turbo Bench Test incorporates the following design:

- ¼ the sample size, increasing exposure to the reaction vessel and rod temperatures.
- Reversed the direction of flow down the Depositor Rod to increase oil contact.
- Used a wire-wound Depositor Rod to control oil exposure time and amount.
- Reduced operating temperatures.
- Removed Nitrous Oxide gas.
- Modified cyclic temperature exposure.
- Redesigned glass mantle for controlling volatile movement within the test.
- Added moist air flow over the Depositor Rod.

TEOST Turbo®: *Test Equipment*

NEW



KING
INFORMATION

TEOST Turbo®: *Test Parameters*

Parameter	33C	MHT®	Turbo®
ASTM Method #	D6335	D7097	D8447
Rod Type	Non-Wound	Wire Wound	Wire Wound
Total Oil Volume	116 mL	~8.5 g	30 mL
Catalyst	6% Fe Naphthenate	MHT®	6% Fe Naphthenate
Gas Induction Location	Reservoir	Upper End Cap	Upper End Cap
Oil Flow	Up the Rod	Down the Rod	Down the Rod
Rod Temp. (Peak)	200°C to 480°C cycle	285° C	290°C w/320°C pk
Reservoir Temp	100°C	N/A	100°C
Test Duration	<2 Hours	24 Hours	18 Hours

TEOST Turbo® Industry Status

ASTM D8447

- ASTM method approved in 2022.


**Industry
Specs.**

- Being considered for industry specification limits.

TEOST® Depositor Rods

33C: unwound

MHT/Turbo: wire-wound





Comments or Questions?

explore our new website
TannasKing.com




EV Lubricant Tests

Conductive Deposit Test
(TanEV CDT)

Wire Corrosion Test
(TanEV WCT)

Brief Introduction

The CDT and WCT data enables a holistic view of the corrosion process and provides the assurance that lubricants which perform well in these tests will be less likely to suffer real world failures in the electrical devices they are supporting.




Coming
Soon

Conductive Deposit Test (CDT)

The formation of a deposit that is electrically conductive is problematic in electric motor designs as this creates a pathway for current to flow away from the intended route and lead to immediate failure of the device.

CDT provides a quantitative measure of conductive deposit or conductive film formation as a function of time and temperature.

Evaluates the potential risk of hardware failures caused by the deposition of conductive corrosion products on electrical components.



Conductive Deposit Test (CDT)



The formation of a deposit that is electrically conductive is problematic in electric motor designs as this creates a pathway for current to flow away from the intended route and lead to immediate failure of the device.

CDT provides a quantitative measure of conductive deposit or conductive film formation as a function of time and temperature.

Evaluates the potential risk of hardware failures caused by the deposition of conductive corrosion products on electrical components.

Specially designed PCB with gaps between copper traces, powered by a 5v DC signal, supplied to the voltage loops. Resistance is tracked during the duration of the test.

PCB lowered into the 20 ml of test fluid at 150C for 500+ hrs. with both liquid and vapor phases being measured for the rate of change in resistance.



Conductive Deposit Test (CDT)

The formation of a deposit that is electrically conductive is problematic in electric motor designs as this creates a pathway for current to flow away from the intended route and lead to immediate failure of the device.


CDT provides a quantitative measure of conductive deposit or conductive film formation as a function of time and temperature.

Evaluates the potential risk of hardware failures caused by the deposition of conductive corrosion products on electrical components.


Specially designed PCB with gaps between copper traces, powered by a 5v DC signal, supplied to the voltage loops. Resistance is tracked during the duration of the test.

PCB lowered into the 20 ml of test fluid at 150C for 500+ hrs. with both liquid and vapor phases being measured for the rate of change in resistance.


Failing Board



Passing Board



Conductive Deposit Test (CDT)



- 4-hour moving average of resistance drops below 50,000 ohms
- Oil fails in Liquid phase at 490 hours
- Oil fails in Vapor phase at 370 hours



Wire Corrosion Test (WCT)

Motor windings / Gaps / Connections and Circuit Boards



- Wide temperature range (80C – 150C)
- Vapor corrosion is a concern


WCT enables a quantitative and mechanistic evaluation of the corrosion tendency of a lubricant in both the solution and vapor phases over a wide range of temperatures (80 to 150C).

Oxidation and the ensuing chemical reactions can cause depletion of critical copper components such as motor windings and interconnects in the EV powertrain.



Wire Corrosion Test (WCT)

400 mL sample,
1 m of Cu wire
(64 μ m diameter)




WCT enables a quantitative and mechanistic evaluation of the corrosion tendency of a lubricant in both the solution and vapor phases over a wide range of temperatures (80 to 150C).

Oxidation and the ensuing chemical reactions can cause depletion of critical copper components such as motor windings and interconnects in the EV powertrain.

Nominal conditions are 130C at 72 hrs. with a constant current applied to the test wire. Test temp and duration can be varied to determine corrosion rate as a function of temperature and time.



Wire Corrosion Test (WCT)



Copper Loss
(Angstroms) vs.
Time



CDT & WCT Industry Status

ASTM

- Both have open working groups in Subcommittee D02.09. Methods expected in late 2023.

Industry

- Written as US Army MILSPEC/FTM and awaiting approval.
- Global EV specifications are being developed via industry OEM groups.
- Included in SAE J3200 for electrified drivetrains.

TANAS

KING
HYDRATION

Tannas Instruments

Precision Laboratory Instruments
for testing Lubricants & Fluids

 <p>TBS Viscometer <i>High-Temperature High-Shear Viscosity</i> ASTM D4683, D6616 CEC L-36 IP370 Tapered Bearing Simulator: Rotational viscometer that measures the high-temperature, high-shear rate viscosity of oils from 40°C to 150°C. Both models available with 40+ position AutoSampler capability. Required test for SAE J300; ILSAC GF Specifications; API SM, SN, SP; ASTM D4485; Chinese GB-11121 & GB-11122.</p>	 <p>Quantum* <i>Non-liquid RPVOT Oxidation Test</i> ASTM D2272, D2112, D4742, D7098 IP229 Oxidation Stability: Dry cylinder, RPVOT (RBOT). Evaluates oxidation resistance over a broad range of oils and lubricants. Widely used for base oil comparisons and condition monitoring for turbine oils. Compact, single-test design without dangerous hot bath oil. Research Accessories: Sample Extraction Device and Sample Temperature Probe</p>	 <p>TanEV CDT <i>Electric Vehicle (EV) Conductive Deposits</i> SAE J3200 ASTM Method Under Development Conductive Deposit Test: Designed to determine the tendency of lubricating fluids to form conductive layer deposits on exposed copper motor windings, connectors, and electrical components at elevated temperatures in an Electrified Vehicle (EV) Drivetrain. A specially modified circuit board (PCB) is lowered into 20 ml of test fluid at 150°C for ~500 hours while supplied with a 5V DC signal.</p>	 <p>TFAB <i>Foam Tendency & Stability</i> ASTM D892, D6082, D1881, D7840 IP146 Tannas Foam Air Bath: A non-liquid bath, the TFAB, tests foam additive effectiveness from 24°C to 150°C. Features a quiet, maintenance-free Direct Drive motor, touchscreen, and a six-position carousel for 1000-mL cylinders. Required test for ILSAC GF Specifications; API SM, SN, SP; ASTM D4485; Chinese GB-11121 & GB-11122.</p>
 <p>Direct Cool® II <i>Low-Temperature Pumpability & Gelation Index</i> ASTM D2983 - Proc. D, D5133, D7110, D8210 The non-liquid, thermoelectric cooling Direct Cool II performs ASTM D2983 and the two SBT® tests, measuring low-temperature pumpability and Gelation Index of fresh, sooted, & highly oxidized oils. Provides continuous rheological data over a broad temperature range (+90°C to -40°C). Required test for ILSAC GF Specifications; API SM, SN, SP; JPI-55-56-99; ASTM D4485; Chinese GB-11121; OEM Specs.</p>	 <p>Noack S2* <i>Evaporation Loss Phosphorus Volatility</i> ASTM D5800 SH/T0059 CEC L-40 Noack S2® Volatility Test: Measures the tendency of base & formulated oils to volatilize in service. Eliminates hazardous Wood's Metal and is tunable to lab environment. Collects volatiles to determine Phosphorus Emission Index (PEI). Required test for ILSAC GF Specifications; API SM, SN, SP; ASTM D4485; Chinese GB-11121 & GB-11122.</p>	 <p>TEOST* <i>High-Temperature Deposit Control</i> ASTM D7097, D6335, D8447 Thermo-oxidation Engine Oil Simulation Test: Measures the high-temperature deposit tendencies of engine oils that form under varying high-temperature conditions (turbo-charger and piston ring areas). For TEOST® 33C, MHT® & TEOST Turbo® tests. Required tests for ILSAC GF Specifications; API SM, SN, SP; ASTM D4485; Chinese GB-11121; OEM factory fill.</p>	 <p>SBT+2 <i>Multi-Purpose Low-Temperature Liquid Bath</i> ASTM D5133, D7110, D2983, D4684, D445, D97, D2500, D2386 Meets numerous low-temperature test methods with easily replaceable Insert Modules for each test method. Designed for use with patented SimAir® Test Cells for low-temperature viscosity of lubricants (ASTM D2983). Ideal for low-temperature work with fresh, sooted, or highly oxidized oils, ATF's, hydraulic fluids and fuels.</p>


TANNAS CO.
4800 James Savage Rd.
Midland, MI 48642 USA
 TannasKing.com
 tannas@savantgroup.com
 +1 989 496 2309
 +1 989 496 3438

 ISO 9001
 © 2023 Tannas Co. | 10/2023 | Printed in USA



King Instruments

*Precision Laboratory Instruments
for testing Lubricants & Fluids*



King MRV TP-1
Mini-Rotary Viscometer
ASTM D3829, D4684, D6821, D6896

Low-Temperature Pumpability & Yield Stress: Determines borderline pumping temperatures of oils and lubricants with patented *direct refrigeration* technology. Features a small bench-top footprint, 10 sample capacity, and removable test cells for ease-of-use.

Required test for SAE J300; ILSAC GF Specifications; API SM, SN, SP; ASTM D4485; JPI 55-42-04; Chinese GB-11121; OEM factory fill.



King BLB
Liquid Bath
ASTM D2983, D97, D2500, D5853 | IP267

Low-Temperature Viscosity: Measures low temperature viscosity of lubricants using the patented SimAir® technique with a small, less costly, *constant* temperature liquid bath. Only technique for independent sample analyses to eliminate batch sample testing. Features a 12-position carousel.

BLB 701 Model: +30°C to -40°C
BLB 702 Model: +30°C to -70°C



King BLB-DIN
BLB-DIN Liquid Bath
DIN 51398 | ASTM D2983 | IP267

Low-Temperature Viscosity: Similar to D2983, the German Standard DIN 51398 determines the apparent viscosity of gear oils and related fluids at low-temperature. Features programmable control, digital display to 0.01°C, and cooling at 1°C/min to -55°C.

The BLB-DIN can be adjusted to run either the D2983 test (using SimAir® Cells) or the DIN standard.



SimAir® Test Cells
*Patented Test Cell
Innovation for ASTM D2983*
ASTM D2983 | IP267

SimAir® Test Cells offer simple, precise, and more efficient data acquisition than Air Baths. Used exclusively in the King BLB and Tannas SBT®+2 liquid baths, the SimAir® Glass Stator design incorporates an insulating chamber between two glass walls to simulate the cooling profile of the original cold air cabinet. Each Test Cell functions with its own independent cooling profile, and can be added or removed from the bath at any time.



KING REFRIGERATION, INC.
4800 James Savage Rd.
Midland, MI 48642 USA

TannasKing.com
tannas@asavantgroup.com

ISO 9001

+1 989 698 5500
+1 989 496 3438



© 2023 King Refrigeration | 03/2023 | Printed in USA



Thank You for your attention!

explore our new website
TannasKing.com

